

Reducing The Time To Diagnosis For Spinal Muscular Atrophy

Zhili Lin, Fen Guo, Yihao Ou, Abhinav Mathur, Zeqiang Ma, and Madhuri Hegde



BACKGROUND

- SMA is the most common neurodegenerative disease in childhood with an incidence of 1 in 6,000 to 1 in 10,000. SMA is caused by deleterious changes in the *SMN1* gene, with a deletion of exon 7 being the most common pathogenic event. Homozygous deletion of exon 7 can be found in approximately 95% of SMA cases, whereas the other 5% are compound heterozygous of this deletion.
- Fragment analysis using the AmpliEx® PCR/CE *SMN1/2* Plus Kit.** The assay is based on PCR and capillary electrophoresis. In addition to *SMN1/2* copy number determination, the assay also detect the presence/absence of gene duplication maskers and modifier c.*3+80T>G, c.*211_*212del, and c.859G>C.
- Identification of SMA positive (*SMN1*=0) and *SMN1* carrier statuses using WGS data.** A bioinformatic workflow was developed and validated for the *SMN1* copy number determination through uniquely mapped reads on exon 7 of *SMN1* gene using the WGS data. Median read depth of *SMN1/SMN2* exon 7 for WGS samples from previous runs are computed and used for normalization.

Validation for the Fragment Analysis Assay

A total of 20 unique samples with known *SMN1/2* copy numbers and /or known variants (15 coriell samples and 5 whole blood clinical samples) selected initially for validation. Two were excluded due to lack of enough gDNA for further testing. All the test results are as expected.


| Results summary of analysis of Coriell and patient samples with known outcomes. | | | | | | | | | |
|---|-------------------------|-------------------------|-------------|----------------|-----------|-------------|--------------|---------------|---------------|
| Identifier | <i>SMN1</i> Copy Number | <i>SMN2</i> Copy Number | c.*3+80 T>G | c.*211_*212del | c.859G >C | Hybrid Peak | As expected? | Analysis Date | Plate Barcode |
| SMAVAL0101 | 0 | 2 | No | No | No | | Yes | 6/23/2020 | C9H11A9H |
| SMAVAL0102 | 0 | 3 | No | No | No | | Yes | 6/23/2020 | C9H11A9H |
| SMAVAL0103 | 1 | ≥4 | No | No | No | | Yes | 6/23/2020 | C9H11A9H |
| SMAVAL0104 | 1 | 1 | No | No | No | | Yes | 6/23/2020 | C9H11A9H |
| SMAVAL0105 | 0 | 3 | No | No | No | | Yes | 6/24/2020 | C9H11A9F |
| SMAVAL0106 | 0 | 2 | No | No | No | | Yes | 6/23/2020 | C9H11A9H |
| SMAVAL0107 | 2 | 0 | No | No | No | | Yes | 6/23/2020 | C9H11A9H |
| SMAVAL0108 | 3 | 0 | Yes | Yes | No | | Yes | 6/25/2020 | C9H11A9E |
| SMAVAL0109 ¹ | ≥4 | 0 | Yes | Yes | No | <i>SMN1</i> | Yes | 6/23/2020 | C9H11A9H |
| SMAVAL0110 | ≥4 | 1 | Yes | Yes | No | | Yes | 6/19/2020 | C9H11A79 |
| SMAVAL0111 | 2 | 2 | No | No | Yes | | Yes | 6/19/2020 | C9H11A79 |
| SMAVAL0112 | 0 | 3 | No | No | No | | Yes | 6/24/2020 | C9H11A9F |
| SMAVAL0113 | 0 | 3 | No | No | No | | Yes | 6/24/2020 | C9H11A9F |
| SMAVAL0114 | 1 | 2 | No | No | No | | Yes | 6/19/2020 | C9H11A79 |
| SMAVAL0115 | 2 | 2 | No | No | No | <i>SMN2</i> | Yes | 6/24/2020 | C9H11A9F |
| SMAVAL0116 | QNS | | | | | | | | |
| SMAVAL0117 | 1 | 3 | No | No | No | | Yes | 6/25/2020 | C9H11A9E |
| SMAVAL0118 | 2 | 2 | No | No | No | | Yes | 6/24/2020 | C9H11A9F |
| SMAVAL0119 | QNS | | | | | | | | |
| SMAVAL0120 | 2 | 0 | No | No | No | | Yes | 6/24/2020 | C9H11A9F |

¹ 2 copies of *SMN1* and 2 copies of *SMN1* hybrid

Identification of SMA (*SMN1*=0) and *SMN1* carrier statuses using WGS data

1) Determination of reference range for *SMN1* copy number

A total of 76 samples were analyzed using the bioinformatic tool, *SMN1* copy numbers were confirmed by MLPA or fragment analysis, and reference ranges were determined for SMA positive (*SMN1*=0), SMA carrier (*SMN1*=1), and no clinical significance (*SMN1*≥2).



2) Identify SMA patients (*SMN1*=0) using WGS data

Three DBS samples were used. The bioinformatic tool identified all 3 as *SMN1*=0.

| Index | ACCESSION | Sex | Bait | Sample type | <i>SMN1</i> _EST_COPY | <i>SMN2</i> _EST_COPY | Orthogonal Method | <i>SMN1</i> _REF_COPY | <i>SMN2</i> _REF_COPY |
|-------|-------------|-----|------|-------------------|-----------------------|-----------------------|-------------------|-----------------------|-----------------------|
| 1 | V2020048296 | F | WGS | Dried Blood Spots | 0 | 2.39 | MLPA | 0 | 2 |
| 2 | V2020074696 | M | WGS | Dried Blood Spots | 0 | 3.11 | MLPA | 0 | 3 |
| 3 | V2020141248 | F | WGS | Dried Blood Spots | 0 | 2.31 | MLPA | 0 | 2 |

3) Identify carrier statuses of *SMN1* using WGS data

WGS samples with estimated *SMN1* copy # within the reference range for SMA carrier were reflexed to the MLPA confirmation test. Here are examples of several clinical samples.

| Sample ID | Estimated <i>SMN1</i> Copy# by WGS | MLPA Confirmed Copy # |
|-----------|------------------------------------|-----------------------|
| WGS #1 | 1.13 | 1 |
| WGS #2 | 1.04 | 1 |
| WGS #3 | 0.83 | 1 |
| WGS #4 | 0.71 | 1 |
| WGS #5 | 1.16 | 2 |
| WGS #6 | 0.86 | 1 |
| WGS #7 | 1.33 | 1 |
| WGS #8 | 1.09 | 1 |
| WGS #9 | 1.33 | 2 |

CONCLUSION

- Currently, this laboratory performs qPCR assay for population-based newborn screening for SMA, bioinformatic analysis of WGS data for identification of SMA patients (*SMN1*=0) and *SMN1* carrier statuses, fragment analysis and MLPA for confirmation and diagnostic testing for determination of *SMN1/2* copy numbers.
- The combination of above assays reduce the time to diagnosis of SMA.